

JAMES A-D Waden



NEU

Länge /D

normal 35–40 cm
lang 41–45 cm

Umfangmasse (cm)

Grösse	cB	cD
XXsmall	18 – 20	23 – 30
XXsmall PLUS	18 – 20	27 – 34
Xsmall	20 – 22	27 – 34
Xsmall PLUS	20 – 22	31 – 38
Small	22 – 24	31 – 38
Small PLUS	22 – 24	34 – 41
Medium	24 – 26	34 – 41
Medium PLUS	24 – 26	37 – 44
Large	26 – 29	37 – 44
Large PLUS	26 – 29	41 – 48
XLarge	29 – 32	41 – 48
XLarge PLUS	29 – 32	42 – 50

Produktlinie

JA2 **JAMES**
ohne Haftrand

Fuss

geschlossen

Farben

black 99, carbon 88,
cyanide 89